



Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India

P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

• Faridabad • Sonapat • Kolkata • Mumbai • Pune • Ahmedabad
• Chennai • Bangalore • Hyderabad

www.prestogroup.com



M 95

Moisture Analyzer



Combining Speed and Precision in One Moisture Analyzer – Test More Samples in Less Time

The M95's halogen heating system distributes heat uniformly over the sample and speeds up the drying process — delivering fast and precise measurements. Designed for tool-free cleaning, the M95's components can be easily removed for quick and easy maintenance. Operating the M95 is intuitive with the icon-driven menu navigation on the touchscreen display.

Standard Features:

- **Increase Efficiency in your Lab with Fast Heating Halogen Technology**

The MB95's precisely controlled halogen heating dries samples quickly and ensures uniform heating to yield fast, precise and repeatable results of 0.01% readability. Four drying profiles and five switch-off criteria help to perform customized sample tests.

- **Define your Drying Parameters and Store Up to 20 Methods for Easy Recall**

MB95 can help you store up to 20 methods and 200 total results, for statistical analysis of measurements.

- **Effortless Maintenance with the Tool-Free Cleaning Design**

Designed with effortless maintenance in mind, the MB95 features a removable pan handler and heating chamber components that make cleaning quick and easy—no tools needed.

- **Touchscreen Icon-Based Menu Navigation Guides Users through Every Step**

Icon-based menu navigation on the easy-to-read touchscreen display makes operating the M95 intuitive. Perform routine sample testing in three easy steps —no training required.