



Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India

P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

• Faridabad • Sonapat • Kolkata • Mumbai • Pune • Ahmedabad
• Chennai • Bangalore • Hyderabad

www.prestogroup.com



GSM Balance Digital



Model No. GSM - 20

Table Top High Accuracy balance

- Capacity : 200gm
- L.C. : 0.01gm
- Pan size : 120 mm
- Display : high contrast LCD display
- Power : AC adaptor (Included) or 4 AA battery (not included)
- Calibration : 200 gm brass weight (Included)
- Linearity : ± 0.01 g
- GSM least count : 1 GSM = (weight x 100)

Optional : Air proof shield cover available on request.

* Other capacity 400gm, 600gm available on request

* Models with higher accuracy also available

Feature :

- Superior shock resistance & overload protection
- Software lock switch & integrated transportation lock
- Adjustable footrest with level indicator.

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-

					
Box Compression Tester	Scuff Resistance Tester	Tensile Tester Digital	Spectrocolorimeter	Cobb Sizing Tester	GSM Round Cutter

www.prestogroup.com

